

# Exhibit A

8

5/22/01

Objective: To coat a monolayer of amino-silane on organic layer on stent to function as a primer prior to polymer coating. This experiment was conducted by Peisen Chang and J. Sunkin.

## Preparation of Amino Silane

0.5ml Amino Silane T2910 (5% Solution)  
9.0 ml Methanol HPLC grade  
0.5 ml Water (De-ionized Sterile)  
10.0 ml Total

Amino Silane Vendor: United Chemical Technology  
Chemical Name: Trimethoxy-silylpropyldiethylamine Triamine  
T2910 Catalog Item. STENTS were immersed for 1 minute in Amino Silane before coating with PCL.

## Preparation of 1% PCL Solution

Weight of PCL = 1.0043 grams  
Volume of THF = 99.0 ml  
Total ~ 100 ml

Lot number of PCL = D99142

## Cleaning of Stents (Total 10 Stents)

1st IPA Sonication = 5 minutes  
1st ~~IPA~~ <sup>DI Water</sup> Sonication = 5 minutes  
2nd IPA Sonication = 5 minutes  
1st DI Water Wash  
2nd DI Water Wash  
NaOH (1N) Sonication = 5 minutes  
DI Water Sonication = 5 minutes  
DI Water Sonication = 5 minutes  
DI Water Sonication = 5 minutes  
Nitrogen Dry (Coat 4 stents during N<sub>2</sub> blow off)  
Heated in Oven (no vacuum) = 15 minutes @ 50°C

R. Sunkin  
5/22/01

5/22/01: Coating of monolayer amino silane as former

Stents were removed from oven after 15 minutes and immersed in 1% PC solution for 10 seconds and placed on Carousel rack for drying. Three coated stents had amino silane treatment (experiment) while the remaining three were control and had no treatment. The Carousel rack was placed in an oven (no vacuum) at 50°C for ~ 1.5 hours. They were cooled for about 5 minutes and inspected under a microscope.

### Microscope observations

It appears that the amino silane coated stents had more wetting of the polymer (PC) than were less were metal. Visible much more coating coverage. The control appeared to be less wettable.

### Coating observations

The 5% solution of amino silane is probably too much for imparting a monolayer. After coating the stent and curing it for 15 minutes we noticed (Peiwen Chang & R. Sundar) that there was white powdery material adhered to stents.

Future concentration of amino silane solution should be targeted for 1% or less. ~~to~~

These experiments are designated as follows:

8-118-17  
-2 } Control  
-3 }

8-118-47  
-5 } with Amino silane T-290  
-6 }

Samples 8-118-1 & 8-118-6 to be submitted to Sam Lab.

R. Sundar  
5/22/01

10

5/25/01

Coating a monolayer of amino silane as primer

SEM observations

SEM observations via several photos show that the amino silane treated stents had less metal exposure on more wetting. However, this could be misleading since the stents were not subjected to uniform dipping by hand. There is some indication that wettability with the use of amino silane is somewhat better but more tests need to be conducted before drawing conclusions.

from  
sundar  
5/25/01

R. Sundar

5/25/01